

VPro™

Quick Start Guide

Your teeth are thanking you - Why VPro?

VPro uses a trusted technology called high frequency vibration (HFV) that stimulates orthodontic tooth movement on a cellular level. When used with your orthodontic treatment for 5 minutes a day, the VPro is clinically proven to:



- ✓ Reduced Treatment Time by up to 64%^{2,4,6}
- ✓ Accelerated Tooth Movement^{2,5,6}
- ✓ Reduced Discomfort^{2,3}
- ✓ Improved Predictability^{2,4}
- ✓ Optimal Aligner Seating¹
- ✓ Assistance in Retainer Wear

From unboxing to smile perfecting - Where to start?

01

Charge
your VPro fully
before your first
session

02

Download
VPro Fastrack™ app
(app not necessary
for device use)

03

Begin
5 minute session
& wear orthodontics
during use

Download the
VPro Fastrack™ App



Light Notifications

Blue (3 blinks) Session Start

Blue (slow blinking) Charging in Progress

Green (1 blink) App has Connected

Green (3 blinks) Session End

Green (slow blinking) Charging Complete

Green (blinks and vibrates) Device Locate

Yellow (2 blinks) Successful App Pairing

Red (slow blinking) Low Battery

Red, Green, Blue (1 blink each) Device Reset

Frequently asked questions

Can I use the VPro more than once a day?

The recommended use is one 5-minute session per day. Please consult with your doctor before exceeding the recommended usage.

How much pressure should I apply?

Bite down gently on the mouthpiece during use. There is no need to bite down harder than as needed to support the device.

How do I clean my mouthpiece?

Remove the mouthpiece from the vibration device and rinse under water, brushing with a toothbrush as necessary.

What if my VPro won't pair or sync with the app?

First, make sure that you have your Bluetooth connection turned on. Next, make sure the VPro is charged and in close proximity to your mobile device. Close and reopen the app. If you are still experiencing difficulty, press and hold the button on the VPro for 10 seconds to reset the device.

